

Scientific contribution from Lonza – Screening and classification of plate structures for enhanced mixing performance

FlowPlate® Reactor Technology to Control Multi-phase Reactions in Continuous-Flow Processing

The control of multi-phase reactions (liquid-liquid (LL) or gas-liquid) is of prime importance to enable broad application of flow technologies. The patented FlowPlate® Reactor Technology answers this challenge via a simple approach based on a deep understanding of reactor design and fluid mechanics. (Figure 1): During the development of a basic reactor structure for multi-phase applications, it was quickly realized that curvature is detrimental to mass transfer and, as a result, needs to be minimized. This work is nicely described in the following publication:

<http://www.sciencedirect.com/science/article/pii/S1385894716305228>

As seen in Figure 2, curvatures generate centrifugal forces which tend to separate phases, not mix them. Despite this, many microreactors or static mixers are designed using curvatures. Clearly this is a pitfall that needs to be avoided when dealing with multi-phase applications.



Figure 2. A typical unwanted Parallel-Flow pattern where the denser phase remains close to the wall

A key objective when developing a multi-phase structure is to generate the so-called drop flow (see Figure 1) at low pressure drop. This is exactly the case with the LL-Rhombus or LL-Triangle structures, though a preference is observed for the LL-Triangle. This is described in “Microreactor mixing-unit design for fast liquid—liquid reactions”

<http://akademai.com/doi/abs/10.1556/1846.2016.00026>.

The next key question is how does such a structure behave with various solvents? This question is answered by Plouffe et al. who describe that interfacial tension has a substantial impact on mass transfer: the higher the value, the lower the mass transfer:

<http://www.sciencedirect.com/science/article/pii/S1385894715014084>

Finally, this research would not be complete without a scale-up study. Scale-up can be ensured at constant energy dissipation rate (pressure drop) when operating in the proper flow regime (e.g. Drop /Dispersed-Flow; see also Figure



Figure 1. FlowPlate patented LL-Triangle-Structure that shows a desirable Drop-Flow

3) The important conclusions are summarized here:

<http://www.sciencedirect.com/science/article/pii/S0009250915007861>




Figure 3. FlowPlate® Lab with LL-Triangle-Structure for small flow rate applications (A– 1–15 mL/min) and a scalable version at ten times this basic flow rate (B– 15-150 mL/min)

A video is also available showing a good-established drop flow:

<https://www.youtube.com/watch?v=5-VEQyAPCBY>

If you cannot access the publications linked to this article, please send us an email at: microreactor@lonza.com

If you have any questions, we will be pleased to answer them by phone, email or in a personal meeting. Visit us under www.ehrfeld.com to obtain an initial impression of our technology.

Or meet us in person at the next event:

9th Symposium on Continuous Flow Reactor Technology for Industrial Applications
14th-16th of November in Barcelona <http://www.flowchemistrytks.com/>

In case of further questions, please do not hesitate to contact us:

info@ehrfeld.com

+49 6734 919300

Best regards,
Ehrfeld Mikrotechnik GmbH

Mikroforum Ring 1, 55234 Wendelsheim, Phone: +49 (0)6734 919 300, info@ehrfeld.com
Geschäftsführung: Dr.-Ing. Joachim Heck, Sitz der Gesellschaft: Wendelsheim, Amtsgericht Mainz HRB 33094
www.ehrfeld.com

If you do not wish to receive e-mails from us please send an e-mail with the subject 'unsubscribe' to anne.kaaden@ehrfeld.com